


8

Retractors






Fig. 2 = 22 X 15 mm
Fig. 1 = 19 X 15 mm




Fig. 1 = 19 X 15 mm
Fig. 2 = 22 X 15 mm

Parker
DC08-001-01 = Set
DC08-001-11 = Fig. 1
DC08-001-12 = Fig. 2
13.5 cm/5 1/4"


Fig. 2 = 27 X 25 mm
Fig. 1 = 24 X 25 mm



Parker
DC08-003-02 = Set
DC08-003-21 = Fig. 1
DC08-003-22 = Fig. 2
18 cm/7"

Fig. 1 = 26 X 25 mm
Fig. 2 = 29 X 25 mm


Fig. 2 = 25 X 24 mm
Fig. 1 = 21 X 24 mm



Parker-Mott
DC08-005-00 = Set
DC08-005-01 = Fig. 1
DC08-005-02 = Fig. 2
16 cm/6 1/4"

Fig. 1 = 68 X 24 mm
Fig. 2 = 72 X 24 mm


34 X 26 mm
52 X 33 mm



Collin-Hartmann
DC08-007-01 = 14 cm/5 1/2"
DC08-007-02 = 16 cm/6 1/4"


Fig. 1 = 24 X 20 mm
Fig. 2 = 29 X 26 mm
Fig. 3 = 35 X 29 mm
Fig. 4 = 27 X 24 mm

Fig. 1 = 30 X 28 mm
Fig. 2 = 36 X 36 mm
Fig. 3 = 38 X 43 mm
Fig. 4 = 42 X 28 mm



Roux
DC08-009-00 = Set Fig. 1-3
DC08-009-01 = Fig. 1 = 14 cm/5 1/2"
DC08-009-02 = Fig. 2 = 15 cm/6"
DC08-009-03 = Fig. 3 = 17 cm/6 3/4"
DC08-009-04 = Fig. 4 = 18 cm/7"

Fig. 1 = 24 X 20 mm
Fig. 2 = 29 X 26 mm
Fig. 3 = 35 X 39 mm
Fig. 4 = 27 X 24 mm



Roux-Magnus
DC08-011-00 = Set Fig. 1-3
DC08-011-01 = Fig. 1 = 14 cm/5 1/2"
DC08-011-02 = Fig. 2 = 15 cm/6"
DC08-011-03 = Fig. 3 = 17 cm/6 3/4"
DC08-011-04 = Fig. 4 = 18 cm/7"

Fig. 1 = 30 X 28 mm
Fig. 2 = 36 X 36 mm
Fig. 3 = 38 X 43 mm
Fig. 4 = 42 X 28 mm






Fig. 2 = 30 X 10 mm
Fig. 1 = 26 X 10 mm




Farabeuf
DC08-013-00 = Set
DC08-013-01 = Fig. 1
DC08-013-02 = Fig. 2
12 cm/4 3/4"


Fig. 2 = 27 X 16 mm
Fig. 1 = 23 X 16 mm



Farabeuf
DC08-015-00 = Set
DC08-015-01 = Fig. 1
DC08-015-02 = Fig. 2
12 cm/4 3/4"




Mini-Roux
DC08-017-06




Baby-Roux
DC08-019-13
13 cm/5 1/4"

Fig. 2 = 47 X 13 mm
Fig. 1 = 42 X 13 mm



Mathieu
DC08-021-00 = Set
DC08-021-01 = Fig. 1
DC08-021-02 = Fig. 2
20 cm/8"



Fig. 2 = 25 X 24 mm
Fig. 1 = 21 X 24 mm




Parker-Langenbeck
DC08-023-00 = Set
DC08-023-01 = Fig. 1
DC08-023-02 = Fig. 2
21 cm/8 1/4"






Goelet
DC08-025-00
19 cm/7 1/2"




Anderson-Adson
DC08-027-20
19.5 cm/8"



Weitlaner
DC08-029-20
20 cm/8"




Weitlaner
DC08-033-10 sharp
DC08-035-10 blunt
10.5cm/4 1/4"




Weitlaner

	sharp	cm	in.	blunt
DC08-037-13		13.0	5	DC08-039-13
DC08-037-16		16.5	6 1/2	DC08-039-16




Weitlaner-Loktite
DC08-041-10 sharp
DC08-043-10 blunt
10.5cm/4 1/4"




Weitlaner-Loktite

	sharp	cm	in.	blunt
DC08-045-13		13.0	5	DC08-047-13
DC08-045-16		16.5	6 1/2	DC08-047-16




Gobble Surgical®



Retractors




Weitlaner-Wullstein
DC08-049-13 sharp
DC08-051-13 blunt
13 cm/5"




Weitlaner
DC08-053-26 sharp
DC08-055-26 blunt
26 cm/10 1/4"




Weitlaner
DC08-057-27 sharp
DC08-059-27 blunt
27 cm/10 1/2"




Milligan
13 cm/5"
DC08-061-13 blunt





Cloward





Cloward








Adson
26.5 cm/10 1/2"



Weitlaner-Baby
DC08-075-13
13 cm/5"



Adson-Baby
DC08-077-16 sharp
DC08-079-16 blunt
16 cm/6 1/2"



Logan
DC08-081-04
4 cm/1 1/2"


Finsen
DC08-082-05 sharp
DC08-083-05 blunt
5 cm/2"





Adson
DC08-084-32
32 cm/12 1/2"




Beckmann-Eaton
DC08-085-32
32 cm/12 1/2"




Beckmann
DC08-086-31
31 cm/12 1/4"




Jansen
DC08-087-10
10 cm/4"
sharp




Jansen
DC08-088-10
10 cm/4"
with 1 sharp and
2 blunt prongs each




Jansen
DC08-089-10
10 cm/4"
blunt




Finsen
DC08-090-07 sharp
DC08-091-07 blunt
7 cm/2 3/4"




Finsen-Mod. Martin
DC08-093-07 sharp
DC08-095-07 blunt
7 cm/2 3/4"



Buckley



Alm





1/2




1/2

cm	in.	Sharp	Sharp	Blunt	Semi-sharp	Blunt
6	2 1/4	DC08-097-06				
7	2 3/4		DC08-101-07	DC08-103-07		
8	3 1/4		DC08-099-08			
10	4		DC08-101-10	DC08-103-10	DC08-105-10	DC08-107-10



Kocher
DC08-109-16
16 cm/6 1/4"

Crotti
DC08-111-18
18.5 cm/7 1/4"



Dartigues
DC08-113-26



Gelpi-Loktite
DC08-115-18
(Perineum)




Gelpi-Loktite
DC08-117-18
(Perineum)





DC08-119-17
17 cm/6 3/4"



Rigby
DC08-121-18
(Rektum, Vagina)




DC08-123-18
18 cm/7"






Mayo-Adams
DC08-125-17




Houzel
DC08-127-16



Baby-Collin
DC08-129-00
complete, with 2 pairs of blades




Collin
DC08-131-22




Collin
DC08-133-00 = 38 X 60 mm/38 X 55 mm
DC08-135-00 = 38 X 80 mm/50 X 75 mm

Measurements
see below



Exchangeable Blades





Collin-Loktite
DC08-137-00




Robin-Masse

	Lateral Blades	Central Blades	Opening width
DC08-139-00	40 X 80 mm	40 X 55 mm	160 mm


Ricard
DC08-141-00 = 80 X 90 mm/80 X 90 mm 300 mm




Gosset
DC08-143-00





Gosset
DC08-145-00




Gosset
DC08-147-00




Gosset
DC08-149-00




Balfour-Baby
DC08-151-00





Balfour
DC08-153-00




Balfour
DC08-155-00



Balfour
DC08-157-00




Balfour
DC08-159-00



Balfour
DC08-161-01 = 45 X 80 mm
DC08-161-02 = 70 X 85 mm
DC08-161-03 = 80 X 90 mm




Doyen
DC08-163-02 = 55 X 095 mm
DC08-163-02 = 55 X 105 mm
DC08-163-03 = 55 X 115 mm
DC08-163-04 = 55 X 125 mm





Qyervain
DC08-165-00




O'Sullivan-O'Connor
DC08-167-00






DC08-169-03
10.5 cm/4 1/4"
10 X 22 mm



DC08-169-06
10.5 cm/4 1/4"
15 X 55 mm




DC08-175-40
40 x 55 mm

DC08-175-50
50 x 65 mm

DC08-173-28

Kirschner
DC08-173-00



Kirschner
DC08-171-00

DC08-175-50
50 x 65 mm

DC08-175-40
40 x 55 mm

DC08-171-24
240 x 240 mm

DC08-175-40
40 x 55 mm



DC08-175-50
50 x 65 mm

DC08-171-00 - Kirschner Abdominal Retractor, complete, consisting of:

- DC08-171-24 - Frame, rigid, 240 x 240 mm
- DC08-175-40 - 2 Hook-on Retractors, 40 x 55 mm
- DC08-175-50 - 2 ditto, 50 x 65 mm
- DC08-171-24 - Frame only, rigid, 240 x 240 mm

DC08-173-00 - Kirschner Abdominal Retractor, folding, large size, complete, consisting of:


- DC08-173-28 - 1 Frame, folding type, large size, 280 x 280 mm
- DC08-175-40 - 2 Hook-on Retractors, 40 x 55 mm
- DC08-175-50 - 2 ditto, 50 x 65 mm
- DC08-173-28 - Frame only, folding type, 280 x 280 mm



Franz (complete)



DC08-179-00
with 4 blades

DC08-181-00
with 8 blades




- DC08-179-00** - **Franz** Abdominal Retractor, folding, complete, consisting of:
- DC08-179-30 - 1 frame, folding, 300 x 300 mm
 - DC08-175-50 - 2 hook-on retractors, 50 x 65 mm
 - DC08-177-38 - 2 exchangeable blades 38 x 6 mm
 - DC08-179-30 - Frame, folding, only, 300 x 300 mm

- DC08-181-00** - **Franz** Abdominal Retractor, folding, complete, consisting of:
- DC08-179-30 - 1 frame, folding, 300 x 300 mm
 - DC08-175-50 - 2 hook-on retractors, 50 x 65 mm
 - DC08-175-60 - 2 ditto, 60 x 85 mm
 - DC08-177-38 - 2 exchangeable blades, 38 x 60 mm
 - DC08-177-45 - 2 ditto, 45 x 80 mm



Temple-Fay
DC08-183-00



Frazier
DC08-185-00



Doyen
DC08-187-25
25 cm/10"



Kocher
DC08-189-22
22 cm/9"







Kocher

	Total length		Blade mm
	cm	in.	
DC08-191-01	25	10	80 X 45
DC08-191-02	25	10	80 X 55
DC08-191-03	25	10	80 X 65



Tuffier
DC08-193-23
23 cm/9"







DC08-203-03
19 cm/7 3/4"




28 cm/11 1/4"




DC08-205-20




DC08-205-25




DC08-205-30





DC08-215-02
75 X 90 mm


1/2



1/2

Davidson
16 cm/6 1/4"
(Scapula)

DC08-215-01
75 X 50 mm



DC08-217-03
82 X 28 mm

DC08-217-02
73 X 28 mm



DC08-217-01
45 X 25 mm


1/2

1/2

1/2



Heaney
27 cm/10 1/2"







Sauerbruch
22.5 cm/9"


Harrington
DC08-227-01 = Fig. 1
DC08-227-02 = Fig. 2
32 cm/12 1/2"
(Splanchnicus)





Brunner
DC08-033-08 - 18
25 cm/10"

1/2



115 x 27mm
115 x 22mm

105 x 27mm

95 x 20mm

$\frac{1}{4}$

$\frac{1}{4}$

$\frac{1}{4}$


$\frac{1}{2}$

Simon
DC08-239-22 = 115 x 22mm
DC08-239-27 = 115 x 27mm
28cm/11"

Martin
DC08-237-27
105 x 27mm
26cm/10 $\frac{1}{4}$ "

Modell Tübingen
DC08-235-95
95 x 20mm
28cm/11"

$\frac{1}{2}$



DC08-241-18 180 x 30mm

DC08-241-16 160 x 30mm

DC08-241-14 140 x 30mm

DC08-241-13 130 x 30mm

DC08-241-12 120 x 25mm

DC08-241-10 100 x 25mm


DC08-241-08 80 x 20mm



$\frac{1}{4}$

$\frac{1}{2}$

Brunner
25cm/10"



$\frac{1}{2}$
Hösel
25cm/10"



DC08-243-62
60 x 20 mm

$\frac{1}{4}$



DC08-243-83
80 x 30 mm

$\frac{1}{4}$



DC08-243-03
100 x 30 mm

$\frac{1}{4}$



DC08-243-23
120 x 30 mm

$\frac{1}{4}$




DC08-243-43
140 x 30 mm

$\frac{1}{4}$



DC08-243-44
140 x 40 mm

$\frac{1}{4}$



38 mm



1/2

Deaver
DC08-245-07
30 cm/12"

50 mm



1/2


Deaver
DC08-245-08
30 cm/12"



75 mm







1/2

Deaver
DC08-245-09
30 cm/12"









14 X 17 mm



DC08-259-01
DC08-261-01
21.5 cm/8 1/2"

20 X 22 mm



DC08-259-02
DC08-261-02
22 cm/8 3/4"


26 X 30 mm



DC08-259-03
DC08-261-03
23.5 cm/9 1/4"



Middeldorpf
DC08-259-01-
DC08-261-03




49 X 38 mm
28 X 20 mm

Fig. 1 ▶

◀ Fig. 2

36 X 28 mm
63 X 49 mm

Richardson-Eastman
DC08-263-00 = Set 1 + 2
DC08-263-01 = Fig. 1 = 26 cm/10 1/4"
DC08-263-02 = Fig. 2 = 28 cm/11"



1/2

Richardson
24 cm/9 1/2"




28 X 20 mm

DC08-265-01




36 X 28 mm

DC08-265-02




44 X 38 mm

DC08-265-03



1/2


Richardson
DC08-267-24
24 cm/9 1/2"



65 X 50 mm

1/2

Kelly
DC08-269-26
26 cm/10 1/4"



15-345-21 = 60 X 20 mm

DC08-271-01 = 35 X 7 mm

DC08-271-02 = 35 X 10 mm

DC08-271-03 = 35 X 15 mm

1/2

Langenbeck
21 cm/8 1/4"





1/2


Kocher
DC08-273-01 = 40 X 10 mm
DC08-273-02 = 40 X 15 mm
DC08-273-03 = 40 X 20 mm
22 cm/8 3/4"



Middeldorpf
DC08-275-01-
DC08-275-03







DC08-287-01 # 1
25 x 6 mm



DC08-287-02 # 2
25 x 8 mm



DC08-287-03 # 3
35 x 8 mm



1/4 DC08-287-04 # 4
35 x 11 mm




DC08-287-05 # 5
35 x 15 mm



DC08-287-06 # 6
40 x 11 mm



DC08-287-07 # 7
55 x 11 mm



DC08-287-08 # 8
70 x 14 mm



DC08-287-09 # 9
80 x 12 mm





1/4 DC08-287-10 # 10
80 x 16 mm



1/2

Kocher-Langenbeck
21.5cm/8 1/2"

also available in special sizes



78 X 40 mm

1/2

DC08-288-02
57 X 25 mm

1/2

DC08-288-01
40 X 18 mm

1/2

DC08-288-03
23 cm/9 1/4"




DC08-289-01
40 X 12 mm 1/2

DC08-289-02
55 X 11 mm 1/2

1/2

Kocher
21 cm/8 1/4"



1/2


DC08-291-01
19 X 12 mm

1/2

DC08-291-02
21 X 14 mm

1/2

Kocher
19 cm/7 1/2"



28 X 6 mm 1/4

1/2


Lahey
DC08-293-19
19 cm/7 1/2"



12 X 14 mm 1/4

1/2

Little
DC08-295-20
20 cm/8"



1/3

Flexible

DC08-297-01 = 33 cm/13" 30 mm
DC08-297-02 = 33 cm/13" 40 mm
DC08-297-03 = 33 cm/13" 50 mm



1/3


Haberer

DC08-299-00 = 20 cm/8" 17/25 mm
DC08-299-01 = 28 cm/11" 25/30 mm
DC08-299-02 = 30 cm/12" 37/45 mm



1/3

DC08-301-12 = 200 X 12mm
DC08-301-17 = 200 X 17mm
DC08-301-25 = 250 X 25mm



1/3

Reverdin

DC08-303-29
29 cm/11 1/2"



1/3


Reverdin

DC08-305-28
28 cm/11"





1/3

DC08-307-28
28 cm/11"
scoop 50X 45 mm




1/3

DC08-309-33
33 cm/13"
85 mm ø
100 ccm



DC08-311-03
12cm/4 3/4"



DC08-311-05
12cm/4 3/4"



O'Conner
DC08-313-12
12cm/4 3/4"



Cushing
DC08-315-14
14cm/5 1/2"



Cloward
DC08-317-16
16.5cm/6 1/2"



Niro
DC08-319-01
17cm/6 3/4"



Niro
DC08-319-02
17cm/6 3/4"




Niro
DC08-319-03
17cm/6 3/4"



Niro
DC08-319-04
17cm/6 3/4"



DC08-321-06
13.5cm/5 1/4"




Miller
13.5cm/6 1/2"



DC08-323-09
sharp




DC08-323-19
blunt



Gobble Surgical®

Fine Wound and Nerve Hooks



Kleinert-Kutz
16cm/6 1/4"

DC08-325-03 DC08-325-05 DC08-325-07



Cushing
DC08-327-14
14.5cm/5 3/4"



Frazier
DC08-329-13
13cm/5"

DC08-329-18
18cm/7"




Frazier
DC08-331-13
13cm/5"

DC08-331-18
18cm/7"



Tonnis
DC08-333-15
15cm/6"



Crile
DC08-335-14
14.5cm/5 3/4"



Crile
DC08-335-20
20cm/8"



Cushing

19cm/7 1/2"
DC08-337-01 DC08-337-02


28cm/11"
DC08-339-01 DC08-339-02



Adson
DC08-341-20
20cm/8"



Adson
DC08-343-20
20cm/8"



Gobble Surgical®

Wound and Trachea Retractors, Mucous Hooks



Jackson
BC08-345-13
13cm/5"



Jackson
BC08-345-15
15cm/6"




News
BC08-347-15
15cm/6"



Cottle
BC08-349-16
16cm/6 1/4"




Freer
BC08-351-15
15cm/6"



Iterson
BC08-353-16
16cm/6 1/4"



Iterson
BC08-355-16
16cm/6 1/4"



Bose
BC08-357-01
16cm/6 1/4"




Bose
BC08-357-02
16cm/6 1/4"



Bose
BC08-357-11
16cm/6 1/4"



Bose
BC08-357-12
16cm/6 1/4"



1/2

Kilner (Barsky)
DC08-359-15
15cm/6"



1/2

Hartmann
DC08-361-16
16cm/6 1/2"



Fig. 1



1/2

Gillies
DC08-363-01
18cm/7"




Fig. 2



1/2

Gillies
DC08-363-02
18cm/7"




1/2

Gillies
DC08-363-52
18cm/7"



1/2

Sachs
DC08-365-21
21cm/8 1/4"



1/2

Schoenborn
DC08-367-21
21cm/8 1/4"



1/4

Lukens
DC08-369-17
17cm/6 3/4"





1/4

Jackson
DC08-371-17
17cm/6 3/4"



1/4

Cope
DC08-373-18
18cm/7"



Day
16cm/6 1/2"





Allport
DC08-377-16
16cm/6 1/2"



Yankauer
DC08-379-15
15cm/6"



Wagner
14cm/5 1/2"



Lucae
14 cm/5 1/2"




Knapp
DC08-385-13
13cm/5"




Rollet
DC08-387-02
12.5 cm/4 3/4"




Rollet (Blair)
DC08-387-04
13.5 cm/5"




Paul
DC08-389-03
12.5 cm/4 3/4"



Bergmann
DC08-391-07
14cm/5 1/2"




Graefe (Tyrell)
12cm/4 3/4"



Guthrie
DC08-397-01
12cm/4 3/4"

Guthrie
DC08-399-02
13cm/5"





Emmet
22cm/8 3/4"



Simon
DC08-403-01
22cm/8 3/4"

Simon
DC08-403-02
22cm/8 3/4"



O'Connor
DC08-405-16
16cm/6 1/4"



Guthrie
16cm/6 1/4"



DC08-407-01



DC08-407-02



Graefe
16 cm/6 1/4"



DC08-409-01
sharp



DC08-409-11
blunt



DC08-411-01



DC08-411-02



DC08-411-03



DC08-411-04



DC08-413-14
14cm/5 1/2"



DC08-415-01




DC08-415-02




DC08-415-03



Cushing
14.5 cm/5 1/2"




Volkmann
23 cm/9"




sharp

1/4



1/2

DC08-417-01



1/2

DC08-417-02



1/2

DC08-417-03




1/2

DC08-417-04



1/2

DC08-417-06



blunt

1/4



1/2

DC08-419-01



1/2

DC08-419-02



1/2

DC08-419-03



1/2

DC08-419-04



1/2

DC08-419-06




1/2

16 cm/6 1/4"



1/4

DC08-421-04




1/4

DC08-423-04



1/2


Cushing
20 cm/8"



8 mm

1/4

DC08-425-08



14 mm

1/4

DC08-425-14



10 mm

1/4

DC08-425-10



16 mm

1/4

DC08-425-16



12 mm

1/4

DC08-425-12



18 mm

1/4

DC08-425-18



1/2

Cushing
24 cm/9 1/2"



8 mm

1/4

DC08-427-08



14 mm

1/4

DC08-427-14



10 mm

1/4

DC08-427-10



16 mm

1/4

DC08-427-16



12 mm

1/4



DC08-427-12



18 mm

1/4

DC08-427-18



Mini-Langenbeck
16 cm/6 1/4"

10 x 6 mm



DC08-429-10

20 x 6 mm



DC08-429-20

17 x 5 mm




DC08-429-17

22 x 8 mm




DC08-429-22

14 x 8 mm



1/2

DC08-431-16
16 cm/6 1/4"



1/2

DC08-433-02
13.5 cm/5 1/4"



1/2

Klapp
17.5 cm/7"



DC08-435-02



DC08-435-03



DC08-437-01











DC08-437-02

toothed







1/2



Strandell-Stille
19 cm/7 1/2"





Micro retractor with
saddle blade, buttoned





O'Connor
DC08-473-16
16 cm/6 1/4"



Graefe
16 cm/6 1/4"



DC08-475-01
sharp



DC08-475-02
blunt



1/4
DC08-477-01



1/4
DC08-477-02



1/4
DC08-477-03



1/4
DC08-477-04

16 cm/6 1/4"



DC08-479-16
16 cm/6 1/4"



Desmarres
16 cm/6 1/4"



DC08-481-08



DC08-481-14



DC08-481-10




DC08-481-16



DC08-481-12




DC08-481-18




DC08-483-12



DC08-483-14




DC08-483-16



DC08-483-18



16 cm/6 1/4"



Guthrie
16 cm/6 1/4"




DC08-485-01



DC08-485-02



16 cm/6 1/4"



DC08-487-01



DC08-487-02



DC08-487-03



DC08-487-04

sharp



16 cm/6 1/4"



DC08-489-01



DC08-489-02




DC08-489-03



DC08-489-04

blunt



DC08-491-02



DC08-493-02



DC08-491-03



DC08-493-03



DC08-491-04



DC08-493-04



DC08-491-06



DC08-493-06

1/2


Volkman
11.5 cm/4 1/2"




Meyerding
18 cm/7"




DC08-495-01




DC08-495-02




DC08-495-03



DC08-495-04



DC08-495-05



DC08-495-06



DC08-495-06

10 x 4 mm








Senn-Green
DC08-497-01 = 15 cm/6"
DC08-497-02 = 15 cm/6"




DC08-497-01



DC08-497-02





Gobble Surgical®


Retractor




$\frac{1}{2}$
Volkmann
21.5 cm/8 $\frac{1}{2}$ "




$\frac{1}{2}$
Kocher
22 cm/8 $\frac{3}{4}$ "







DC08-525-07
27 cm/10 3/4"




DC08-527-04
40 X 40 mm



DC08-529-05
45 X 50 mm



DC08-531-06
50 X 60 mm



DC08-533-07
70 X 70 mm



Israel
25 cm/10"

