





16

Kidney, Liver, Gall, Urology







1
DC16-001-01
23cm/9"




2
DC16-001-02
22cm/8¾"




3
DC16-001-03
22cm/8¾"




4
DC16-001-04
22cm/8¾"




Randall



5
DC16-001-05
19cm/7"



6
DC16-001-06
19cm/7"





1/4

DeBakey
DC16-003-20
20cm/8"


Cooley
DC16-005-18
18.5cm/7 1/4"

Cooley
DC16-005-19
19.5cm/7 3/4"


Renal Artery Clamps- also used for subclavin and carotid artery




Stille
DC16-007-22
22cm/8 ¾"




Stille
DC16-009-24
24cm/9 ½"



Wertheim
DC16-011-25
25cm/10"



Herrick
DC16-013-23
23cm/9"




Pean



DC16-015-20
20cm/8"

DC16-015-22
22cm/8 ¾"


DC16-015-24
24cm/9 ½"




Guyon
DC16-017-23
23cm/9"





DC16-019-20





Blake
DC16-021-20
20 cm / 8"





Mayo-Blake
DC16-023-20
20 cm / 8"





Mayo-Blake
DC16-025-20
20 cm / 8"




Czerny
DC16-027-25
24.5 cm / 9 3/4"



Turner-Warwick
DC16-029-18
18 cm / 7"



Turner-Warwick
DC16-031-19
19 cm / 7 1/2"




Gobble Surgical®




Gall Stone Forceps, Gall Duct- and Cystic Forceps




Mixer
DC16-033-22
22 cm / 8 3/4"




Desjardins
23 cm / 9"




Desjardins
DC16-039-21
21cm/8 1/4"



Nissen
DC16-041-22
22.5cm/9"




Green
DC16-043-20
20cm/8"



Gray

DC16-045-22
22cm/8 3/4"


DC16-045-23
23cm/9"



Moynihn


DC16-047-23
23cm/9"

DC16-047-24
24cm/9 1/2"





Gobble Surgical®


Gall Stone Probes, Gall Duct Dilators



- DC16-049-38 - Gall Duct Probe, with stilete, malleable
- DC16-053-28 - DC16-057-28 - with eccentric olives
- DC16-059-28 - DC16-063-28 - Gall Stone Probes, elastic
- DC16-065-01 - DC16-065-13 - Gall Duct Dilators, Copper shaft, malleable











Mod. Muenchen
30 cm / 12"
elastic




Desjardins
30 cm / 12"
elastic




Luer-Koerte
32 cm / 12 3/4"
malleable





Gobble Surgical®

Gall Stone Scoops



1/4

Fergusson
DC16-087-21
21.5 cm / 8 1/2"




1/4

Fergusson
DC16-087-22
22 cm / 8 3/4"



1/4

Fergusson
DC16-087-23
23.5 cm / 9 1/4"




1/2

Mayo
DC16-089-27
28 cm / 11"





1/2

DC16-091-34
34 cm / 12 1/2"




1/2


Moore
DC16-093-28
28 cm / 11"



DC16-095-16
16 cm / 6 1/4"
Spermatic Cord Hook, blunt



Millin
DC16-097-23
23 cm / 9"
Capsule Forceps,
for retropubic prostatectomy




1/2

Millin
DC16-099-23
23 cm / 9"
prostatic Lobe Holding Forceps,
for anterior lobes, 9:9 teeth




Millin
DC16-101-24
24 cm / 9 1/2"
Ligature Carrying Forceps




Gobble Surgical®


Prostatic Retractors, Circumcision Instrument




Young
DC16-103-21
21 cm / 8 1/4"




Young
DC16-105-22
22 cm / 8 3/4"



Young
DC16-107-21
21 cm / 8 1/4"
for bulb




Young
DC16-109-22
22 cm / 8 3/4"




Circumcision Instrument
DC16-111-08 - DC16-111-32

- DC16-111-08 - 8 mm ø, baby size
- DC16-111-11 - 11 mm ø, child size
- DC16-111-13 - 13 mm ø, child size
- DC16-111-16 - 16 mm ø, adolescent size
- DC16-111-21 - 21 mm ø, adolescent size
- DC16-111-26 - 26 mm ø, men's size
- DC16-111-32 - 32 mm ø, men's size



Gobble Surgical®

Dilating Bougies, Set of Dilating Sounds



Dittel
6-10 mm
18-30 Charriere

Dittel
Charriere (French)
No. 8-30

Van Buren
Charriere (French)
No. 8-36



Guyon
Charriere
No. 10-30

Guyon
Charriere
No. 10-30


Filiform
31-610-12

Set	DC16-115-00	DC16-117-00	DC16-119-00	DC16-121-00
single	DC16-113-00	DC16-115-08 - 31-680-30	DC16-117-08 - 31-699-36	DC16-119-10 - DC16-119-30

	Scale mm Ø	I	II	Dittel	Van Buren curved	Guyon without int. thread	Guyon with int. thread
Scale I Charriere (French) No. 8-36 rising by 1/3 mm per number	2 2/3	8		DC16-115-08	DC16-117-08		
	3	9		DC16-115-09			
	3 1/3	10	20	DC16-115-10	DC16-117-10	DC16-119-10	DC16-121-10
	3 2/3	11	22	DC16-115-11		DC16-119-11	DC16-121-11
	4	12	24	DC16-115-12	DC16-117-12	DC16-119-12	DC16-121-12
	4	13	26	DC16-115-13		DC16-119-13	DC16-121-13
	4	14	28	DC16-115-14	DC16-117-14	DC16-119-14	DC16-121-14
	5	15	30	DC16-115-15		DC16-119-15	DC16-121-15
	5	16	32	DC16-115-16	DC16-117-16	DC16-119-16	DC16-121-16
	5	17	34	DC16-115-17		DC16-119-17	DC16-121-17
	6	18	36	DC16-115-18	DC16-117-18	DC16-119-18	DC16-121-18
	6	19	38	DC16-115-19		DC16-119-19	DC16-121-19
Scale II Benique No. 20-26, rising by 1/3 mm per two numbers	6	20	40	DC16-115-20	DC16-117-20	DC16-119-20	DC16-121-20
	7	21	42	DC16-115-21		DC16-119-21	DC16-121-21
	7	22	44	DC16-115-22	DC16-117-22	DC16-119-22	DC16-121-22
	7	23	46	DC16-115-23		DC16-119-23	DC16-121-23
	8	24	48	DC16-115-24	DC16-117-24	DC16-119-24	DC16-121-24
	8	25	50	DC16-115-25		DC16-119-25	DC16-121-25
	8	26	52	DC16-115-26	DC16-117-26	DC16-119-26	DC16-121-26
	9	27	54	DC16-115-27		DC16-119-27	DC16-121-27
	9	28	56	DC16-115-28	DC16-117-28	DC16-119-28	DC16-121-28
	9	29	58	DC16-115-29		DC16-119-29	DC16-121-29
	10	30	60	DC16-115-30	DC16-117-30	DC16-119-30	DC16-121-30
	10	32			DC16-117-32		
	10	34			DC16-117-34		
	12	36			DC16-117-36		



DC16-123-00
male and female



DC16-125-00
male, female and
for children

Table of measurements

Scale mm ø	Charriere (French)
2	6
2 1/3	7
2 2/3	8
3	9
3 1/3	10
3 2/3	11
4	12
4 1/3	13
4 2/3	14
5	15
5 1/3	16
5 2/3	17
6	18
6 1/3	19
6 2/3	20
7	21
7 1/3	22
7 2/3	23
8	24


DC16-127-06	DC16-129-06
DC16-127-07	DC16-129-07
DC16-127-08	DC16-129-08
DC16-127-09	DC16-129-09
DC16-127-10	DC16-129-10
DC16-127-12	DC16-129-12
DC16-127-13	DC16-129-13
DC16-127-14	DC16-129-14
DC16-127-15	DC16-129-15
DC16-127-16	DC16-129-16
DC16-127-17	DC16-129-17
DC16-127-18	DC16-129-18
DC16-127-19	DC16-129-19
DC16-127-20	DC16-129-20
DC16-127-21	DC16-129-21
DC16-127-22	DC16-129-22
DC16-127-23	DC16-129-23
DC16-127-24	DC16-129-24



DC16-127-06 - DC16-127-24
Charriere (French) No. 6-24
27 cm / 10 1/2"
Male Catheters




DC16-129-06 - DC16-129-24
Charriere (French) No. 6-24
15 cm / 6"
Female Catheters Combined,
for pocket sets



Gobble Surgical®

Foreign Body Forceps, Meatotome
Secretion Scoops, penis Seizing Forceps, Meatus Clamps



Lengths of shanks

Mathieu

DC16-131-20 = 20 cm/8"
DC16-131-28 = 28 cm/11"

Urethral Foreign Body Forceps, serrated jaws


Mathieu

DC16-133-20 = 20 cm/8"
DC16-133-28 = 28 cm/11"

Urethral Foreign Body Forceps, serrated jaws
sharp, punch jaws

Sachse

DC16-135-12
Meatotome




Guttman

DC16-137-26
Secretion Scoops flexible,
may be held in the open flame,
without holes
DC16-139-26
with holes
26.5 cm / 10 1/2"




Jusbasic

DC16-141-19
19 cm / 7 1/2"
Penis Seizing Forceps




Strauss

DC16-143-11
Meatus Clamp, rubber covered



Judd-Masson
DC16-145-00

1/2



Judd-Masson
DC16-147-00

1/2



Thomson-Walker
DC16-149-00

1/2



1/2



Legueu
DC16-151-00

1/2

- DC16-147-00 - for supra pubic prostatectomy, with 2 swivel blades
- DC16-149-00 - self-retaining, with 4 detachable blades, center spatula and retractor
- DC16-151-00 - with 2 detachable blades and center spatula



35 mm

1/2

Legueu
DC16-153-27
27 cm / 10 1/2"



56 mm

1/2



Kocher
DC16-155-26
26 cm / 10 1/4"





67 mm

1/2

Kocher
DC16-157-26
26 cm / 10 1/4"





Millin
DC16-161-00
Meatus Clamp
(Bladder Neck Spreader)

Gosset
DC16-163-00
Bladder Retractor